LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

INSTRUCTION FOR USE-GEBRAUCHSANWEISUNG

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Lens Delivery System is a combination of Cartridge, Cushion, and Injector. Lens Delivery System is a single-use and sterile device intended to insert one-piece foldable intraocular lens.

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II. DEVICE DETAILS: - GERÄTEDETAILS:

There are various models of lens delivery system varying in their cartridge size 2.2mm to 2.4mm. There are different vendors for lens delivery system. All our IOL's are compatible for all our LDS so depending on client's requirement the LDS are provided.

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There are models with similar specification but different brands this is due to marketing purpose

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LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

Device Sub- category- Geräte- Unterkategori e	Model Number Modellnummer	Recommended Compatible LDS model-Empfohlenes kompatibles LDS-Modell	
		Diopter Range(-) 5 to (+)24.5D Diopter-Bereich (-) 5 bis (+) 24,5 D	Diopter Range(+) 25 to (+)45D Diopter-Bereich (+) 25 bis (+) 45 D
	YSQFL600ASP*	AQ-S-B-CON22, AQ-S-B-CON24, AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-Flyglide, IA-S-B-MD3, IA-S-B-MD4 , IA-S-B-RT3, IA-S-B-RT4	AQ-S-B-CON24, AQ-S-B-MJ24 IA-S-B-MD4 IA-S-B-RT4
	YSQQ600ASP*/HQL1110606EYA	AQ-S-B-CON22, AQ-S-B-CON24, AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-Flyglide, IA-S-B-MD3, IA-S-B-MD4 , IA-S-B-RT3, IA-S-B-RT4	AQ-S-B-CON24, AQ-S-B-MJ24 IA-S-B-MD4 IA-S-B-RT4
	SQQ600ASP*	IA-S-B-MD4 , IA-S-B-RT3, IA-S-B-RT4	
	SQQ600MCR*+/PBF37UVASP*+/AS68/AS64	AQ-S-B-CON22, AQ-S-B-CON24, AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-Flyglide, IA-S-B-MD3, IA-S-B-MD4 , IA-S-B-RT3, IA-S-B-RT4	AQ-S-B-CON24, AQ-S-B-MJ24 IA-S-B-MD4 IA-S-B-RT4
	YSQQ600MCR*+/AS68Y		
	YSQFL600ASP/HFL2125605NYA	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide

LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

Device Sub- category- Geräte- Unterkategori e	Model Number Modellnummer	Recommended Compatible LDS model-Empfohlenes kompatibles LDS-Modell	
		Diopter Range(-) 5 to (+)24.5D Diopter-Bereich (-) 5 bis (+) 24,5 D	Diopter Range(+) 25 to (+)45D Diopter-Bereich (+) 25 bis (+) 45 D
	YSQQ600ASP/ AY4/LDS AY4/HPL1110605CRA	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide
	SQQ600ASP/HQL1110606CRA		
	PBF37UVASP	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide
	UVF600-130Q	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide

LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

Device Sub- category- Geräte- Unterkategori e	Model Number Modellnummer	Recommended Compatible LDS model-Empfohlenes kompatibles LDS-Modell	
		Diopter Range(-) 5 to (+)24.5D Diopter-Bereich (-) 5 bis (+) 24,5 D	Diopter Range(+) 25 to (+)45D Diopter-Bereich (+) 25 bis (+) 45 D
	CBFY33SUVASP+/AS66Y/HLL1125605NYA	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide
	SQA600ASP/CBF33SUVASP+/AS66	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide
Hydrophilic Foldable Multifocal IOL	YSQFL600DF/HFL2125605NYD	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide

LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

Device Sub-category- Geräte- Unterkategorie	Model Number Modellnummer	Recommended Compatible LDS model-Empfohlenes kompatibles LDS-Modell	
		Diopter Range(-) 5 to (+)24.5D Diopter-Bereich (-) 5 bis (+) 24,5 D	Diopter Range(+) 25 to (+)45D Diopter-Bereich (+) 25 bis (+) 45 D
	YSQQ600DF/M AY4	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide
	PBFY37MF/HQL1110605NYD	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide
Hydrophilic Foldable Toric IOL	CBFY33UVT#	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide
	PBFY37UVQT#/TAY4	AQ-S-B-MJ22, AQ-S-B-CON22, IA-S-B-MD3, IA-S-B-RT3, AQ-S-B-Flyglide	AQ-S-B-MJ22, AQ-S-B-MJ24, AQ-S-B-CON22, AQ-S-B-CON24, IA-S-B-MD3, IA-S-B-MD4, IA-S-B-RT3, IA-S-B-RT4, AQ-S-B-Flyglide

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The lens delivery system is one of the most effective and useful methods for inserting intraocular lens (IOL) during cataract surgery. In comparison to manual insertion of the lens, IOL implanted using an injector has demonstrated small incision width and higher wound stability¹. With the advancements in intraocular lens materials and designs, the concept of dispensing an IOL as pre-packaged, ready-to-insert IOL has been widely adopted by manufacturers². The primary function of the injector system is to deliver the lens seamlessly into a capsular bag with minimal inflammation and a smaller incision.

A variety of options are available for the user to choose the right injector system such as semi-preloaded or fully pre-loaded delivery systems. This technical documentation consists of a lens delivery system that is supplied with a sterile injector and cartridge set. Parts of the lens delivery system are purchased from different suppliers and are assembled, packaged, and sterilized by Omni Lens Pvt. Ltd.

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Cartridge, made of polypropylene, is a device used to hold the lens in place and penetrate the eye for lens implantation. It is available in different sizes with a varying range of tip diameters. The design of the cartridge comprises flaps, guide rail, and the tip.

The function of the guide rail is to hold the lens in place, in a foldable manner if possible. The flaps, also known as wings, ensure the complete fitting of the lens inside the cartridge. The tip penetrates inside the eye for transient time and helps deliver the lens in the capsular bag.

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LENS DELIVERY SYSTEM-LINSENABGABESYSTEM


Figure 1.1 Cartridge-Kartusche

Injector primarily consists of the injector body, plunger, and spring as shown in figure 1.2 below. The injector body is an outer body made of ABS (TR) material. And Plunger is made from Poly Carbonate. The function of the plunger is to push the lens forward with the help of the spring attached to it. A flange is designed on the injector body to hold the injector while pushing the plunger. The spring attached to the plunger tip absorbs extra resistance from the force applied through the plunger to enable smooth movement of the lens.

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The cushion is a soft silicone rubber that prevents the lens from being damaged when pushed through the plunger. The cushion is placed on the tip of the plunger. Omni Lens Pvt. Ltd. carries out assembling of the cushion into the injector followed by packaging and sterilization of injector and cartridge.

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
a. Polypropylene-Polypropylen

Polypropylene (PP), also known as polypropene, is a thermoplastic polymer used in a wide variety of applications. It is produced via chain-growth polymerization from the monomer propylene. Polypropylene belongs to the group of polyolefins and is partially crystalline and non-polar. Its properties are similar to polyethylene, but it is slightly harder and more heat resistant. It is a white, mechanically rugged material and has a high chemical resistance. Bio-PP is the bio-based counterpart of polypropylene (PP).

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LENS DELIVERY SYSTEM-LINSENABGABESYSTEM


Polypropylene is in many aspects similar to polyethylene, especially in solution behaviour and electrical properties. The methyl group improves mechanical properties and thermal resistance, although the chemical resistance decreases. The properties of polypropylene depend on the molecular weight and molecular weight distribution, crystallinity, type and proportion of comonomer (if used) and the isotacticity. In isotactic polypropylene, for example, the methyl groups are oriented on one side of the carbon backbone. This arrangement creates a greater degree of crystallinity and results in a stiffer material that is more resistant to creep than both atactic polypropylene and polyethylene.

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b. Silicone- Silikon

A silicone or polysiloxane is a polymer made up of siloxane ($-\text{R}_2\text{Si}-\text{O}-\text{SiR}_2-$, where R = organic group). They are typically colorless oils or rubber-like substances. Silicones are used in sealants, adhesives, lubricants, medicine, cooking utensils, thermal insulation, and electrical insulation. Some common forms include silicone oil, silicone grease, silicone rubber, silicone resin, and silicone caulk.

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LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

More precisely called polymerized siloxanes or polysiloxanes, silicones consist of an inorganic silicon–oxygen backbone chain ($\text{---Si-O-Si-O-Si-O---}$) with two organic groups attached to each silicon center. Commonly, the organic groups are methyl. The materials can be cyclic or polymeric. By varying the ---Si-O--- chain lengths, side groups, and crosslinking, silicones can be synthesized with a wide variety of properties and compositions. They can vary in consistency from liquid to gel to rubber to hard plastic. The most common siloxane is linear polydimethylsiloxane (PDMS), a silicone oil. The second-largest group of silicone materials is based on silicone resins, which are formed by branched and cage-like oligosiloxane

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
c. Polycarbonate-Polycarbonate

Polycarbonates (PC) are a group of thermoplastic polymers containing carbonate groups in their chemical structures. Polycarbonates used in engineering are strong, tough materials, and some grades are optically transparent. They are easily worked, molded, and thermoformed. Because of these properties, polycarbonates find many applications. Polycarbonates do not have a unique resin identification code (RIC) and are identified as "Other", 7 on the RIC list. Products made from polycarbonate can contain the precursor monomer bisphenol A (BPA).

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Polycarbonate is a durable material. Although it has high impact-resistance, it has low scratch-resistance. Therefore, a hard coating is applied to polycarbonate eyewear lenses and polycarbonate exterior automotive components. The characteristics of polycarbonate compare to those of polymethyl methacrylate (PMMA, acrylic), but polycarbonate is stronger and will hold up longer to extreme temperature. Thermally processed material is usually totally amorphous, and as a result is highly transparent to visible light, with better light transmission than many kinds of glass.

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d. Stainless Steel-Edelstahl

Stainless steel, originally called rustless steel, is any one of a group of ferrous alloys that contain a minimum of approximately 11% chromium, a composition that prevents the iron from rusting and also provides heat-resistant properties. Different types of stainless steel include the elements carbon, nitrogen, aluminium, silicon, sulfur, titanium, nickel, copper, selenium, niobium, and molybdenum. Specific types of stainless steel are often designated by their AISI three-digit number, e.g., 304 stainless.

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All patients who need an implantation of a foldable intraocular lens.

- Vision loss
- Decrease in quality of life
- Depression
- Low visual acuity
- Difficulty in performing visual task

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VII. MODE OF ACTION: WIRKUNGSWEISE:

Injector/Cartridge (Lens delivery system) is used to deliver the intraocular lens in capsular bag during the cataract surgery.

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LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

X. METHOD OF STERILIZATION: STERILISATIONSMETHODE:

Lens Delivery Systems are supplied dry, in a package, terminally sterilized with Ethylene oxide and must be opened under aseptic conditions.

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Reducing, as far as possible, the risk of use error due to the ergonomic features of the device and the environment in which the device is intended to be used (design for patient safety).

1. Avoidance of potential IOL loading errors
2. Reduced operation time
3. Provided better operating conditions
4. Provide effective and safe IOL delivery
5. provide a clinical improvement in IOL delivery for patients undergoing cataract surgery
6. Gives safe, effective and predictable deliver IOLs in the eye
7. Reducing the potential for lens contamination and postsurgical complications

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XII. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

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Link for the availability to SSCP: Link to be provided after notified body acceptance of SSCP.

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Store & transport between 5°C to 40°C, Keep away from sun light.

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The use of a lens delivery system is essential for the implantation of intraocular lens. It consists of cartridge, injector and cushion.

Lens delivery system is supplied in single/Regular Pack.

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
01/08/2025


LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

4. Open the blister carefully & take out the cartridge from its blister pack in a sterile environment.

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LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

7. Place the foldable IOL in such a manner that the lens is centred and front haptic pointing left.

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As with any surgical procedure, there is risk involved. The following non-exhaustive list specifies the complications that have been associated with the implantation of IOLs:

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II. RESIDUAL RISKS- RESTRISIKEN

The finished device is having the Residual Risks such as

LENS DELIVERY SYSTEM-LINSENABGABESYSTEM

- Delivery Failure
- Allergic Reaction

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XXII. EXPIRATION DATE INFORMATION: INFORMATIONEN ZUM VERFALLSDATUM:

Sterility is guaranteed unless the pouch is damaged or opened. The expiration date is clearly indicated.

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XXIII. RETURN GOODS POLICY: RÜCKGABEPOLITIK:

Omni Lens Pvt. Ltd. accepts returned LDS for exchanges only in case of manufacturing defect. No cash refunds will be issued. To return LDS, you must first obtain a Return authorization number from customer services department. No returned goods will be accepted without proper authorization number. Returned LDS should be shipped by traceable method. No credit will be given to lost or damaged LDS shipment. LDS will be replaced as long as they are returned within six months of their original invoice date.

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- The manufacturer guarantees that this product was produced with appropriate care and shall assume no responsibility for incidental or consequential damages, losses or costs that should result directly or indirectly from the use of this product.
- Liability is solely limited to claim-related repairs that must be performed on the product which are clearly not attributed to incorrect handling or the use of lenses not validated for this injector model.
- **Der Hersteller garantiert, dass dieses Produkt mit angemessener Sorgfalt hergestellt wurde und übernimmt keine Verantwortung für mittelbare oder Folgeschäden, Verluste oder Kosten, die direkt oder indirekt auf die Nutzung dieses Produkts zurückzuführen sind.**
- **Die Haftung beschränkt sich ausschließlich auf Gewährleistungsreparaturen, die an dem Produkt vorgenommen werden müssen und eindeutig nicht auf eine falsche Handhabung oder den Einsatz von Linsen zurückzuführen sind, die für dieses Injektormodell nicht validiert wurden.**

Omni Lens Pvt. Ltd.

EU REP

Obelis S.A.
Boulevard General Wahis 53
1030, Brussels, Belgium

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